


Subdermal needle electrodes EMG

Fast, simple and safe in application!!

Needle electrodes for single use

Subdermal needle electrodes present the simplest yet most reliable method of recording EMG signals. Motor nerves (such as the facial nerve) can therefore be safely monitored during surgical interventions via the respective target muscles.



Art. no.: 50-5510

Set of 3 needle electrodes, 20 mm

Single-channel EMG

Packed under sterile conditions P/10

Art. no.: 50-5512

Set of 3 needle electrodes, 13 mm

Single-channel EMG

Packed under sterile conditions P/5